
ACAMO



PIP Endoprosthesis PRESSFIT

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Development and Production
of Ceramic Implants made of
Zirconium Ceramic with a
Bio Active Ceramic Coating

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

- Quick pain relief
- Restoration of good strength
- Maintenance and restoration of mobility
- Very good biocompatibility
- No abrasion
- Cement-free
- Simple OP-Technique

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The operative surgery takes place under plexus anaesthetic, the upper arm is exsanguinated. Long, longitudinal skin incision of approx. 3-4 cm in the middle of the proximal interphalangeal (PIP) joint. Longitudinal opening of the dorsal aponeurosis and the extensor tendon proximal interphalangeal joint, exposure of the base of the middle phalanx.

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Exposure of the joint surfaces. In cases of severe deformation of the base of the middle phalanx peel away the middle phalanx after longitudinally severing it at the base and moving it to the sides. Protection of the ligament base radially and ulnarly at the head of the proximal phalanx.

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Positioning of the sawing gauge proximally at the shaft of the proximal phalanx and resection of the joint surface at the relevant angle of inclination towards palmar. Distal resection of joint surface at the base of the middle phalanx. In the case of deformed joints this is often only possible with the resection of the base of the middle phalanx and the base of the lateral ligament.

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Altogether approx. 8 mm lengths are needed for the two prosthesis parts. Check with the distance jig.

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Position the K-wire into the shaft of the proximal phalanx using image intensification (using a drilling machine or hand-drill chuck) in longitudinal axis. K-Wire 1.0 mm

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The distance and drilling jig with a central guide hole can be used for the exact positioning of the K-wire. Overdrill the K-wire and drill open the shaft of the proximal phalanx with the drill (depending on the chosen shaft thickness with diameter 2.5-2.8 or 3.1)

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The spongiosa is then compacted with the conical compactor (according to the prosthesis size with 2.7 or 3.0). There is a blunt-edged and a sharp-edged rasp available. Rasps/compactors proximally (proximal phalanx).

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In the same way, now the shaft of the middle phalanx is prepared. For this the rasps/milling compactors for 2.5-2.7 and 3.0 shafts are used. Re-check the position and tension by placing the space holders.

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Rasps/milling compactors for the distal part of the joint (middle phalanx).

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The proximal part of the prosthesis is first placed into the shaft of the proximal phalanx with the help of an impactor (concave). Axis and rotation can be controlled with the toggle.

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The distal part of the prosthesis is then put in place after opening up the joint. It is put in place with the convex impactor and controlled by the toggle.

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Check mobility and positioning of the prosthesis using an X-ray image intensifier.

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Reconstruction and refixation of the lateral ligaments and middle phalanxes.

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After the operation place a palmar palm-finger splint at the inclination of the MCP joint of about 60°. Extension of middle and distal phalanx.

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Important information about our orthopaedic implants:

Endoprotheses are intended for the replacement of normal body structures. When using implants, the surgeon should be aware of the following.

1. Selecting the correct implant is extremely important.

The potential for successful surgery is greatly increased by selecting the proper type and size of implant. The size and shape of the human bone place restrictions on the size and shape of the implant, thus limiting its load-bearing capabilities. No implant can be expected to withstand the unsupported stress of full body weight. Endoprotheses are especially dependent upon proper placement and adequate bone support and should not be expected to withstand more than normal stresses and functional body weight.

2. Correct handling of the implant is extremely important.

During the operation the implant should under no circumstances be subjected to blows from hard objects. These and similar treatment can produce defects in the surface finish causing internal stress concentrations which may weaken the implant and lead to possible failure. Only the appropriate instruments should be used during the installation of the implant.

3. No surgical implant should ever be re-used or re-sterilised locally.

Even if a previously used implant appears to be undamaged, it may already have developed internal fatigue which can lead to failure. We urge you to use only new implants of the latest design and construction.

4. Post-operative care is important.

The patient must be made aware of the limitations of his new implant. They must be cautioned that the new implant can only be expected to withstand limited stress until the supporting bones are completely healed. The load-bearing capability of implants cannot be compared to that of a healthy bone.

MOJE – KERAMIK Petersberg

Detailed information regarding important criteria can be supplied upon request. This catalogue was developed for the sale and promotion of our products described and illustrated herein. Reproduction of its contents, or any part thereof, is expressly prohibited. In the event of misuse, we reserve the right to reclaim both catalogue and price list as well as the right to recourse.

All instruments, unless otherwise noted, are made of stainless steel.

The implants are supplied in sterile packing.

Design, Patents, Production, Certification

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articotec GmbH (Langenhagen)
Moje Keramik Implantate GmbH & Co.KG
articotec GmbH